

ELISpot substrate: BCIP/NBT-plus for ALP

Product code:	3650-10
Contents:	Contains BCIP 0.692 mM and NBT 0.734 mM in a 2-amino-2-methyl-1-propanol buffer at pH 9.8. Also contains non-toxic stabilizers.
Applications:	Substrate for alkaline phosphatase conjugates in ELISpot assay.
Instructions:	Filter the ready-to-use substrate solution (BCIP/NBT-plus) through a 0.45 µm filter and add to ELISpot plates after incubation with alkaline phosphatase and subsequent washing. Develop until spots emerge. Stop the colour development by washing extensively in tap water.
Storage:	Shipped at ambient temperature. Store at 4-8°C upon receipt. Adjust to room temperature before use.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com
Updated on 2022-11-22

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.