

Goat anti-mouse IgG+IgM-Alkaline Phosphatase

Product code:	3310-5
Antibody:	Antiserum is raised in goat using purified antigen as immunogen. The antibody is specific for mouse IgG and IgM H+L-chains and was purified using immunoaffinity chromatography. No reactivity against human, horse and bovine serum proteins or mouse non-immunoglobulin serum proteins. Conjugate may cross-react with immunoglobulins from other species.
Enzyme:	Highly purified enzyme from calf intestine with high activity especially suited for enzyme immunoassay purposes.
Conjugation:	Alkaline Phosphatase is conjugated to the H ₂ N-ligand by employing a modified protocol which uses glutardialdehyde as the crosslinking agent. The preparation has been thoroughly depleted of free ligand and enzyme. The conjugate has high residual enzyme activity and antigen affinity.
Quantity:	1 ml
Applications:	Enzyme immunoassays, Western blot techniques.
Working dilution:	Supplied in 0.1 M Tris, 1 mM MgCl ₂ , 0.1 mM ZnCl ₂ , 1% BSA, 0.002 % Kathon CG, pH 7.4. Dilute 1/1000 as a general guideline. Dilution factor may vary with differing assay conditions.
Storage:	Shipped at ambient temperature. On arrival store product at 4-8°C.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Intertek



Intertek

Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the above standards.

MABTECH AB
Box 1233
SE-131 28 Nacka Strand
Sweden
Tel +46 8 716 27 00
Fax +46 8 716 27 01
mabtech@mabtech.com

MABTECH Inc.
M.E.B. Suite 220, 3814 West Street
OH 45227 Cincinnati
USA
Tel: +1 513 871 4500
Fax: +1 513 871 7353
mabtech.usa@mabtech.com
Toll Free: +1 866 ELI-SPOT

MABTECH Australia Pty Ltd
Australia
Tel: +61 3 9470 4704
Fax: +61 3 8678 3216
mabtech.au@mabtech.com

MABTECH AB
Büro Deutschland
Germany
Tel: +49 40 4135 7935
Fax: +49 40 4135 7945
mabtech.de@mabtech.com

MABTECH AB
Bureau de liaison
France
Tel: +33 4 92 38 80 70
Fax: +33 4 92 38 80 71
mabtech.fr@mabtech.com