## Mouse IFN-γ FluoroSpot<sup>BASIC</sup> (550)

PRODUCT CODE: FS1-41-550

Capture antibody: mAb AN18, 1 mg/ml (600 µl)

**Detection antibody:** mAb R4-6A2, biotinylated, 1 mg/ml (50 μl)

Fluorophore conjugate: SA-550 (240 µl)

**Fluorescence enhancer:** Ready-to-use Fluorescence enhancer-II (25 ml)

Supplied in: Monoclonal antibodies are supplied in sterile filtered (0.2 µm) PBS with 0.02% sodium

azide. Streptavidin-conjugate SA-550 and the Fluorescence enhancer contain 0.002%

Kathon CG.

Quantity: Reagents for 4 plates

**Storage:** Shipped at ambient temperature. On arrival, store product refrigerated at 4-8°C.

**Specificity:** Mouse IFN-γ

**Applications:** For enumeration of IFN-γ secreting cells using FluoroSpot. Please find the recommended

assay conditions below.

**Analysis:** A FluoroSpot reader equipped with a filter for Cy3 (excitation 550 nm/emission 570 nm)

is recommended for analysis.

**Protocols:** General protocol for the FluoroSpot assay is provided in the package and is also available

at www.mabtech.com\fluorospot-protocol.

**Assay conditions:** Reagents should not be diluted in advance

Capture antibody	AN18	Dilute to 15 µg/ml
Detection antibody	R4-6A2, biotinylated	Dilute to 1 μg/ml
Fluorophore conjugate	SA-550	Dilute 1:200
Fluorescence enhancer	Fluorescence enhancer-II	50 μl/well
Cell incubation time	12-48 hours in the plate	
Recommended positive	ConA (1-10 μg/ml)	
controls		



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.





Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the above standards.

Mabtech AB (Head Office) Sweden Tel +46 8 716 27 00 mabtech@mabtech.com Augustendalstorget 9 SE - 131 52 Nacka Strand Mabtech Inc. USA Tel: +1 513 871-4500 mabtech.usa@mabtech.com

Mabtech AB, Büro Deutschland Germany Tel: +49 40 4135 7935 mabtech.de@mabtech.com Mabtech AB, Bureau de liaison France Tel: +33 (0)4 92 38 80 70 mabtech.fr@mabtech.com Mabtech Australia Pty Ltd Australia Tel: +61 3 9470 4704 mabtech.au@mabtech.com