

Monoclonal antibody to Hamster Interferon- γ

Antibody:	MTH21
Product code:	3102-3-250
Quantity:	250 μ g
Immunogen:	Recombinant hamster (<i>Mesocricetus auratus</i> golden hamster) IFN- γ
Isotype:	Mouse IgG1
Specificity:	Native and recombinant hamster IFN- γ
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IFN- γ in solution e.g. cell culture supernatants using ELISA and for enumeration of IFN- γ producing cells using ELISpot. MTH21 is recommended as coating mAb in combination with detection mAb MTH29 (product code 3102-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com