

# Monoclonal antibody to Bovine Interleukin-2

<b>Antibody:</b>	MT3B3
<b>Product code:</b>	3111-5-250
<b>Quantity:</b>	250 µg
<b>Immunogen:</b>	Recombinant bovine IL-2
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Native and recombinant bovine IL-2. Cross-reacts with IL-2 from sheep.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For detection of IL-2 in cow and sheep using Flow cytometry.
<b>Protocols:</b>	For further information, please contact Mabtech.

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com