

# Monoclonal antibody to Equine Interferon- $\gamma$

<b>Antibody:</b>	MT166
<b>Product code:</b>	3117-3-1000
<b>Quantity:</b>	1000 $\mu$ g
<b>Immunogen:</b>	Recombinant equine IFN- $\gamma$
<b>Isotype:</b>	Mouse IgG2a
<b>Specificity:</b>	Native and recombinant IFN- $\gamma$ from horse. Cross-reacts with native and recombinant IFN- $\gamma$ from dog and native IFN- $\gamma$ from rhinoceros.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 $\mu$ m) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of equine IFN- $\gamma$ in solution e.g. cell culture supernatants using ELISA and for enumeration of equine IFN- $\gamma$ producing cells using ELISpot. MT166 is recommended as coating mAb in combination with detection mAb MT13 (product code 3117-6).
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com