

Monoclonal antibodies to Porcine Interferon- γ

Antibody:	pIFN γ -I
Product code:	3130-3-1000
Quantity:	1000 μ g
Immunogen:	Recombinant porcine IFN- γ
Isotype:	Mouse IgG1
Specificity:	Native porcine IFN- γ . Recombinant porcine IFN- γ produced in <i>E.coli</i> and human cells. No reactivity with recombinant porcine IFN- γ produced in yeast.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IFN- γ in solution e.g. cell culture supernatants using ELISA and for enumeration of IFN- γ producing cells using ELISpot. pIFN γ -I is recommended as coating mAb in combination with detection mAb P2C11-biotin (product code 3130-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com