

# Monoclonal antibody to Porcine TNF- $\alpha$

<b>Antibody:</b>	MT1P14
<b>Product code:</b>	3135-1-250
<b>Quantity:</b>	250 $\mu$ g
<b>Immunogen:</b>	Recombinant porcine TNF- $\alpha$
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Native and recombinant porcine TNF- $\alpha$
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 $\mu$ m) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of porcine TNF- $\alpha$ in solution e.g. cell culture supernatants using ELISA. MT1P14 is recommended as coating mAb in combination with detection mAb MT859 (product code 3135-6).
<b>Protocols:</b>	Please contact Mabtech for ELISA protocol.

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com