

# ELISpot conjugate: Streptavidin-ALP

<b>Product code:</b>	3310-10-1000
<b>Streptavidin:</b>	Recombinant streptavidin with high activity, especially suited for enzyme immunoassay purposes.
<b>Enzyme:</b>	Highly purified enzyme from calf intestine with high activity especially suited for enzyme immunoassay purposes.
<b>Conjugation:</b>	The preparation has been thoroughly depleted of free ligand and free enzyme. The conjugate has high residual enzyme activity and biotin affinity.
<b>Quantity:</b>	1 ml
<b>Applications:</b>	Optimised for ELISpot. For other enzyme immunoassays please inquire Mabtech.
<b>Working dilution:</b>	Supplied in 0.1 M Tris, 1 mM MgCl <sub>2</sub> , 0.1 mM ZnCl <sub>2</sub> , 1% BSA, 0.002 % Kathon CG, pH 7.4. Dilute 1/1000 as a general guideline. Dilution factor may vary with differing assay conditions.
<b>Storage:</b>	Shipped at ambient temperature. On arrival store product at 4-8°C. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on 2024-02-09

[www.mabtech.com](http://www.mabtech.com)