

ELISA/ELISpot conjugate: Streptavidin-HRP

Product code:	3310-9-1000
Streptavidin:	Recombinant streptavidin with high activity, especially suited for enzyme immunoassay purposes.
Enzyme:	Highly purified enzyme from horseradish with high activity especially suited for enzyme immunoassay purposes.
Conjugation:	The enzyme and streptavidin have been covalently conjugated. The conjugate has high residual enzyme activity and biotin affinity.
Quantity:	1 ml
Applications:	Enzyme immunoassays such as ELISA, ELISpot and Western Blot.
Working dilution:	Supplied in PBS, 1% BSA, 0.002% Kathon CG, pH 7.4 Dilute 1/1000 as a general guideline. Please note that when used for ELISpot the dilution factor may vary significantly depending on the substrate used. <i>HRP-conjugates should not be used with buffers containing sodium azide as this compound will inhibit enzyme activity.</i>
Storage:	Shipped at ambient temperature. On arrival store product at 4-8°C. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2024-02-09

www.mabtech.com