

Monoclonal Antibody to Mouse Interleukin-6

BIOTINYLATED

Antibody:	32C11
Product code:	3361-6-1000
Quantity:	1000 µg
Immunogen:	Recombinant mouse IL-6
Isotype:	Rat IgG2a
Specificity:	Native and recombinant mouse IL-6
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified by affinity chromatography
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of mouse IL-6 in solution, e.g. cell culture supernatants, using ELISA and for enumeration of IL-6 producing mouse cells using ELISpot. 32C11 is recommended as detection mAb in combination with coating mAb MT26B2 (product code 3361-3).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note: for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com