

Monoclonal antibody to Monkey/Human Interferon- γ

Antibody:	MT126L
Product code:	3421M-3-1000
Quantity:	1000 μ g
Immunogen:	Recombinant human IFN- γ
Isotype:	Rat IgG2a
Specificity:	Native and recombinant human IFN- γ . Native IFN- γ from rhesus, cynomolgus and pigtail macaques, baboon, sooty mangabey, African green monkey, aotus, common marmoset and squirrel monkey. Please visit www.mabtech.com for information on reactivity.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IFN- γ in solution e.g. cell culture supernatants using ELISA and for enumeration of IFN- γ producing cells using ELISpot. MT126L is recommended as coating mAb in combination with detection mAb 7-B6-1 (product code 3420-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com