

# Monoclonal antibody to Monkey Interleukin-10

<b>Antibody:</b>	MT3515
<b>Product code:</b>	3430M-3-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Recombinant common marmoset IL-10
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Native and recombinant IL-10 from common marmoset and native IL-10 from rhesus macaque and cynomolgus macaque. Please visit <a href="http://www.mabtech.com">www.mabtech.com</a> for information on reactivity with other monkey species. The antibody cross-reacts with IL-10 from human and horse.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of monkey IL-10 and horse IL-10 in solution e.g. cell culture supernatants using ELISA and for enumeration of monkey IL-10 and horse IL-10 producing cells using ELISpot. MT3515 is recommended as coating mAb in combination with detection mAb MT32AD (product code 3430M-6).
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
[mabtech@mabtech.com](mailto:mabtech@mabtech.com)

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
[mabtech.usa@mabtech.com](mailto:mabtech.usa@mabtech.com)