

# Monoclonal Antibody to Human Interleukin-2

BIOTINYLATED

<b>Antibody:</b>	MT8G10
<b>Product code:</b>	3445-6-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Recombinant human IL-2
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Native and recombinant human IL-2. Native and recombinant IL-2 from cynomolgus macaque. Please visit <a href="http://www.mabtech.com">www.mabtech.com</a> for information on reactivity with other monkey species.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of IL-2 in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-2 producing cells using ELISpot. MT8G10 is recommended as detection mAb in combination with coating mAbs MT2A91/2C95 (product code 3445-3). MT8G10 is also recommended for use in flow cytometry.
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.

# MABTECH

Note: for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
[mabtech@mabtech.com](mailto:mabtech@mabtech.com)

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
[mabtech.usa@mabtech.com](mailto:mabtech.usa@mabtech.com)

Updated on, 2022-11-14

[www.mabtech.com](http://www.mabtech.com)