

Monoclonal antibody to Human Granzyme A

Antibody:	GA29
Product code:	3484-3-1000
Quantity:	1000 µg
Immunogen:	Human Granzyme A
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human Granzyme A.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Affinity chromatography
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of Granzyme A in solution e.g. cell culture supernatants using ELISA. GA29 is recommended as coating mAb in combination with detection mAb GA28 (product code 3484-6).
Protocols:	Please contact Mabtech for ELISA protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com