

# Monoclonal Antibody to Human Interleukin-29 (IFN- $\lambda$ 1)

BIOTINYLATED

<b>Antibody:</b>	MT6G4
<b>Product code:</b>	3570-6-1000
<b>Quantity:</b>	1000 $\mu$ g
<b>Immunogen:</b>	Recombinant human IL-29
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Native and recombinant human IL-29. Cross-reacts with human IL-28A and IL-28B.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 $\mu$ m) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of human IL-29 in solution, e.g. cell culture supernatants, using ELISA. MT6G4 is recommended as detection mAb in combination with IL-29 specific coating mAb MT2H7 (product code 3570-3).
<b>Protocols:</b>	Please contact Mabtech for ELISA protocols.

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on, 2022-11-18

[www.mabtech.com](http://www.mabtech.com)