

Monoclonal antibody to Human CD3

Antibody:	CD3-2
Product code:	3605-1-1000
Quantity:	1000 µg
Immunogen:	Stimulated human leukocytes
Isotype:	Mouse IgG2a
Specificity:	Recognizes the ε chain of human CD3
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS. Contains no preservatives.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	This mAb is useful for activation of T cells. It can be administered in solution or a solid phase. Recommended concentration for optimal cytokine production is 100 ng/ml. It is also useful for studies of CD3 expression in cells and tissue sections. For flow cytometry staining of paraformaldehyde fixed cells, fixation is recommended after staining.
Protocols:	For further information, please contact Mabtech.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com