

FluoroSpot enhancer

Product code:	3641-F10
Contents:	Ready-to-use solution with 0.002% Kathon CG as preservative.
Applications:	The FluoroSpot enhancer will increase spot quality in FluoroSpot assays. Please note that the FluoroSpot enhancer is optimized for use with FluoroSpot assays from Mabtech.
Storage:	Shipped at ambient temperature. Store at 4-8°C upon receipt. 120 ml

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com
Updated on 2022-11-22

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.