

Monoclonal antibody to Monkey/Human Apolipoprotein E

Antibody:	E981
Product code:	3712M-3-1000
Quantity:	1000 µg
Immunogen:	Purified human apolipoprotein E
Isotype:	Mouse IgG2a
Specificity:	Native and recombinant human apolipoprotein E (isoforms apoE2, apoE3 and apoE4). Native apoE from rhesus and cynomolgus macaques. Please visit www.mabtech.com for information on reactivity with other monkey species. The presence of low amounts of detergent is necessary for binding. Tween 20, Triton X100 or NP40 can be used at a concentration of 0.01-0.5%. Do not vortex with detergent added.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of apoE in serum/plasma samples and cell culture supernatants using ELISA. E981 is recommended as coating mAb in combination with detection mAb E887 (product code 3712-6). The antibody is also suitable for immunoprecipitation and Western blot.
Protocols:	Please contact Mabtech for ELISA protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com