

Monoclonal antibody to Human Apolipoprotein M

Antibody:	M5
Product code:	3716-3-250
Quantity:	250 µg
Immunogen:	Recombinant human apoM
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human apoM
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of apoM in solution e.g. serum/plasma samples and cell culture supernatants using ELISA. M5 is recommended as coating mAb in combination with biotinylated detection mAb M7 (product code 3716-6). The antibody is also suitable for immunoprecipitation and Western blot.
Protocols:	Please contact Mabtech for ELISA protocols.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com