

Monoclonal Antibody to Mouse Apolipoprotein A1

BIOTINYLATED

Antibody:	mHDL36
Product code:	3750-6-1000
Quantity:	1000 µg
Immunogen:	Purified mouse apoA1
Isotype:	Rat IgG2a
Specificity:	Mouse apoA1
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of apoA1 in solution e.g. serum/plasma samples and cell culture supernatants using ELISA. mHDL36 is recommended as detection mAb in combination with coating mAb mHDL93 (product code 3750-3). The antibody is also suitable for Western blot.
Protocols:	Please contact Mabtech for ELISA protocols.

MABTECH

Note: for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com