

Monoclonal antibodies to Human IgE

Antibody:	107/182/101 The reagent consists of a combination of three monoclonal antibodies specific for human IgE.
Product code:	3810-7-1000
Quantity:	1000 µg
Immunogen:	Human IgE
Isotype:	Mouse IgG1/IgG2b/IgG2a
Specificity:	Human IgE
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For detection of IgE in e.g. antigen-specific ELISA.
Protocols:	Please contact Mabtech for more information.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2022-11-14

www.mabtech.com