

Monoclonal Antibody to Human IgM

BIOTINYLATED

Antibody:	MT22
Product code:	3880-6-1000
Quantity:	1000 µg
Immunogen:	Human IgM
Isotype:	Mouse IgG2b
Specificity:	Human IgM. Cross-reacts with IgM from rhesus and cynomolgus macaques. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For enumeration of IgM producing cells in ELISpot and for quantification of antigen-specific IgM using e.g. indirect ELISA. In ELISpot mAb MT22 is recommended as detection mAbs in combination with coating mAbs MT11/12 (product code 3880-3).
Protocols:	Please contact Mabtech for ELISpot protocols.

MABTECH

Note: for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com