

# Monoclonal Antibody to Human IgM

BIOTINYLATED

<b>Antibody:</b>	MT22
<b>Product code:</b>	3880-6-250
<b>Quantity:</b>	250 µg
<b>Immunogen:</b>	Human IgM
<b>Isotype:</b>	Mouse IgG2b
<b>Specificity:</b>	Human IgM. Cross-reacts with IgM from rhesus and cynomolgus macaques. Please visit <a href="http://www.mabtech.com">www.mabtech.com</a> for information on reactivity with other monkey species.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For enumeration of IgM producing cells in ELISpot and for quantification of antigen-specific IgM using e.g. indirect ELISA. In ELISpot mAb MT22 is recommended as detection mAbs in combination with coating mAbs MT11/12 (product code 3880-3).
<b>Protocols:</b>	Please contact Mabtech for ELISpot protocols.

# MABTECH

Note: for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
[mabtech@mabtech.com](mailto:mabtech@mabtech.com)

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
[mabtech.usa@mabtech.com](mailto:mabtech.usa@mabtech.com)

Updated on, 2022-11-14

[www.mabtech.com](http://www.mabtech.com)