

Monoclonal antibody to Woodchuck Interferon- γ

Antibody:	MT247
Product code:	3101-3-1000
Quantity:	1000 μ g
Immunogen:	Recombinant woodchuck IFN- γ
Isotype:	Mouse IgG1
Specificity:	Native and recombinant woodchuck IFN- γ
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of woodchuck IFN- γ in solution using ELISA and for enumeration of woodchuck IFN- γ producing cells using ELISpot. MT247 is recommended as coating mAb in combination with detection mAb MT225 (product code 3101-6).
Protocols:	Please contact Mabtech for protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2021-03-18



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com