

Monoclonal Antibody to Woodchuck Interferon- γ

BIOTINYLATED

Antibody:	MT225
Product code:	3101-6-1000
Quantity:	1000 μ g
Immunogen:	Recombinant woodchuck IFN- γ
Isotype:	Mouse IgG1
Specificity:	Native and recombinant woodchuck IFN- γ
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of woodchuck IFN- γ in solution using ELISA and for enumeration of woodchuck IFN- γ producing cells using ELISpot. MT225 is recommended as detection mAb in combination with coating mAb MT247 (product code 3101-3).
Protocols:	Please contact Mabtech for protocols.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tél: +46 8 716 27 00
mabtech@mabtech.com

Updated on, 2021-03-18

Mabtech, Inc.
USA
Tél: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com