

# Monoclonal antibody to Hamster Interferon- $\gamma$

<b>Antibody:</b>	MTH21
<b>Product code:</b>	3102-3-250
<b>Quantity:</b>	250 $\mu$ g
<b>Immunogen:</b>	Recombinant hamster ( <i>Mesocricetus auratus</i> golden hamster) IFN- $\gamma$
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Native and recombinant hamster IFN- $\gamma$
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 $\mu$ m) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For enumeration of IFN- $\gamma$ producing cells using ELISpot. MTH21 is recommended as coating mAb in combination with detection mAb MTH29 (product code 3102-6).
<b>Protocols:</b>	Please contact Mabtech for ELISpot protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on 2021-01-14



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)