

Monoclonal antibody to Bovine Interleukin-2

Antibody:	MT3B3
Product code:	3111-5-250
Quantity:	250 µg
Immunogen:	Recombinant bovine IL-2
Isotype:	Mouse IgG1
Specificity:	Native and recombinant bovine IL-2. Cross-reacts with IL-2 from sheep.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For detection of IL-2 in cow and sheep using Flow cytometry.
Protocols:	For further information, please contact Mabtech.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com