

Monoclonal antibodies to Bovine Interleukin-8 (CXCL8)

Antibody:	MT8H6/8F19 The reagent consists of a combination of two different monoclonal antibodies specific for human IL-8.
Product code:	3114-1-1000
Quantity:	1000 µg
Immunogen:	Recombinant human IL-8
Isotype:	Mouse IgG1
Specificity:	Native and recombinant bovine IL-8. Native and recombinant human IL-8. Cross-reacts with IL-8 from cynomolgus macaque and IL-8 from dog.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For enumeration of IL-8 producing cells using ELISpot. MT8H6/8F19 are recommended as coating mAbs in combination with detection mAb 26E5 (product code 3114-6).
Protocols:	Please contact Mabtech for ELISpot protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com