

Monoclonal antibody to Cotton rat Interferon- γ

Antibody:	MT413
Product code:	3121-3-1000
Quantity:	1000 μ g
Immunogen:	Recombinant cotton rat IFN- γ
Isotype:	Mouse IgG1
Specificity:	Native and recombinant cotton rat IFN- γ
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of cotton rat IFN- γ in solution e.g. cell culture supernatants using ELISA and for enumeration of cotton rat IFN- γ producing cells using ELISpot. MT413 is recommended as coating mAb in combination with detection mAb MT417 (product code 3121-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com