

# Monoclonal antibody to Chicken Interferon- $\gamma$

<b>Antibody:</b>	MT6C2
<b>Product code:</b>	3125-3-250
<b>Quantity:</b>	250 $\mu$ g
<b>Immunogen:</b>	Recombinant chicken IFN- $\gamma$
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Native and recombinant chicken IFN- $\gamma$ .
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 $\mu$ m) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of chicken IFN- $\gamma$ in solution e.g. cell culture supernatants using ELISA and for enumeration of chicken IFN- $\gamma$ producing cells using ELISpot. MT6C2 is recommended as coating mAb in combination with detection mAb MT7C10 (product code 3125-6).
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on 2020-04-14



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)