

Monoclonal Antibody to Chicken Interferon- γ

BIOTINYLATED

Antibody:	MT7C10
Product code:	3125-6-1000
Quantity:	1000 μ g
Immunogen:	Recombinant chicken IFN- γ
Isotype:	Mouse IgG1
Specificity:	Native and recombinant chicken IFN- γ .
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of chicken IFN- γ in solution e.g. cell culture supernatants using ELISA and for enumeration of chicken IFN- γ producing cells using ELISpot. MT7C10-biotin is recommended as detection mAb in combination with coating mAb MT6C2 (product code 3125-3).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tél: +46 8 716 27 00
mabtech@mabtech.com

Updated on, 2020-04-15

Mabtech, Inc.
USA
Tél: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com