

Monoclonal antibodies to Porcine Interleukin-2

Antibody:	MT264/267
Product code:	3140-3-1000
Quantity:	1000 µg
Immunogen:	Recombinant porcine IL-2
Isotype:	Mouse IgG1
Specificity:	Native and recombinant porcine IL-2
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IL-2 in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-2 producing cells using ELISpot. MT264/267 are recommended as coating mAbs in combination with detection mAb MT265 (product code 3140-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-06-22



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com