

Monoclonal antibody to Bovine IgG

Antibody:	MT134
Product code:	3150-3-1000
Quantity:	1000 µg
Immunogen:	Bovine IgG
Isotype:	Rat IgG1
Specificity:	Bovine IgG Cross-reacts with IgG from human (86% cross-reactivity), horse (6% cross-reactivity) and sheep (1% cross-reactivity).
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For enumeration of bovine IgG producing cells using ELISpot. MT134 is recommended as coating mAb in combination with detection mAb MT391 (product code 3150-6). MT134 can also be used as capture mAb in ELISA in combination with detection mAb MT391.
Protocols:	Please contact Mabtech for ELISpot protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com