

# Monoclonal antibody to Mouse Interferon- $\gamma$

|                       |   |
|-----------------------|---|
| <b>Antibody:</b>      | AN18  |
| <b>Product code:</b>  | 3321-3-250  |
| <b>Quantity:</b>      | 250 $\mu$ g   |
| <b>Immunogen:</b>     | Partially purified mouse IFN- $\gamma$  |
| <b>Isotype:</b>       | Rat IgG1  |
| <b>Specificity:</b>   | Native and recombinant mouse IFN- $\gamma$  |
| <b>Concentration:</b> | Supplied at 1 mg/ml in sterile filtered (0.2 $\mu$ m) PBS with 0.02% sodium azide.  |
| <b>Purification:</b>  | Purified from <i>in vitro</i> cultures by protein G affinity chromatography.  |
| <b>Storage:</b>       | Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.   |
| <b>Applications:</b>  | For quantification of IFN- $\gamma$ in solution e.g. cell culture supernatants using ELISA and for enumeration of IFN- $\gamma$ producing cells using ELISpot. AN18 is recommended as coating mAb in combination with detection mAb R4-6A2 (product code 3321-6). |
| <b>Protocols:</b>     | Please contact Mabtech for ELISpot and ELISA protocols.   |



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)