

Monoclonal antibodies to Mouse Interferon- α

Antibody:	MT24A/44A/104
Product code:	3326-1-250
Quantity:	250 μ g
Immunogen:	Recombinant mouse IFN- α
Isotype:	Rat IgG1
Specificity:	Native and recombinant mouse IFN- α . MT24A/44A/104 in combination with antibodies MT9L/14A/113 detect subtypes 1, 2, 4, 5, 6T, 7/10 (129/Sv), 8/6, 9, 11, 12, 13, 14, 15 (A), and B.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IFN- α in solution e.g. cell culture supernatants using ELISA. MT24A/44A/104 are recommended as coating mAbs in combination with detection mAbs MT9L/14A/113 (product code 3326-14).
Protocols:	Please contact Mabtech for ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2021-03-18



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com