

Monoclonal Antibody to Mouse Interleukin-22

BIOTINYLATED

Antibody:	MT231
Product code:	3376-6-1000
Quantity:	1000 µg
Immunogen:	Recombinant mouse IL-22
Isotype:	Mouse IgG1
Specificity:	Native and recombinant mouse IL-22. Cross-reacts with rat IL-22.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For enumeration of IL-22 producing cells using ELISpot. MT231 is recommended as detection mAb in combination with coating mAb MT230 (product code 3376-3).
Protocols:	Please contact Mabtech for ELISpot protocols.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tél: +46 8 716 27 00
mabtech@mabtech.com

Updated on, 2020-03-26

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com