

# Monoclonal antibody to Mouse Interleukin-5

<b>Antibody:</b>	TRFK5
<b>Product code:</b>	3391-3-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Mouse IL-5
<b>Isotype:</b>	Rat IgG1
<b>Specificity:</b>	Native and recombinant mouse IL-5. Native and recombinant human IL-5
<b>Concentration:</b>	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of mouse IL-5 in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-5 producing cells using ELISpot. TRFK5 is recommended as coating mAb in combination with detection mAb TRFK4 (product code 3391-6).
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on 2020-04-07



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)