

Monoclonal antibody to Monkey/Human Interferon- γ

Antibody:	GZ-4
Product code:	3420M-3-1000
Quantity:	1000 μ g
Immunogen:	Recombinant human IFN- γ
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human IFN- γ . Native IFN- γ from rhesus, cynomolgus and pigtail macaques. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Applications:	For quantification of IFN- γ in solution e.g. cell culture supernatants using ELISA and for enumeration of IFN- γ producing cells using ELISpot. GZ-4 is recommended as coating mAb in combination with detection mAb 7-B6-1 (product code 3420-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Intertek



Intertek

MABTECH AB
Box 1233
SE-131 28 Nacka Strand
Sweden
Tel +46 8 716 27 00
Fax +46 8 716 27 01
mabtech@mabtech.com

MABTECH Inc.
M.E.B. Suite 220, 3814 West Street
OH 45227 Cincinnati
USA
Tel: +1 513 871 4500
Fax: +1 513 871 7353
mabtech.usa@mabtech.com
Toll Free: +1 866 ELI-SPOT

MABTECH Australia Pty Ltd
Australia
Tel: +61 3 9470 4704
Fax: +61 3 8678 3216
mabtech.au@mabtech.com

MABTECH AB
Büro Deutschland
Germany
Tel: +49 40 4135 7935
Fax: +49 40 4135 7945
mabtech.de@mabtech.com

MABTECH AB
Bureau de liaison
France
Tel: +33 4 92 38 80 70
Fax: +33 4 92 38 80 71
mabtech.fr@mabtech.com