

# Monoclonal antibody to Monkey/Human Interferon- $\gamma$

<b>Antibody:</b>	MT126L
<b>Product code:</b>	3421M-3-250
<b>Quantity:</b>	250 $\mu$ g
<b>Immunogen:</b>	Recombinant human IFN- $\gamma$
<b>Isotype:</b>	Rat IgG2a
<b>Specificity:</b>	Native and recombinant human IFN- $\gamma$ . Native IFN- $\gamma$ from rhesus, cynomolgus and pigtail macaques, baboon, sooty mangabey, African green monkey, aotus, common marmoset and squirrel monkey. Please visit <a href="http://www.mabtech.com">www.mabtech.com</a> for information on reactivity.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 $\mu$ m) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of IFN- $\gamma$ in solution e.g. cell culture supernatants using ELISA and for enumeration of IFN- $\gamma$ producing cells using ELISpot. MT126L is recommended as coating mAb in combination with detection mAb 7-B6-1 (product code 3420-6).
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
[mabtech@mabtech.com](mailto:mabtech@mabtech.com)

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
[mabtech.usa@mabtech.com](mailto:mabtech.usa@mabtech.com)



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Updated on 2020-03-26

[www.mabtech.com](http://www.mabtech.com)