

Monoclonal antibody to Human Interferon- α (subtype 2)

Antibody:	MT1
Product code:	3423-3-1000
Quantity:	1000 μ g
Immunogen:	Human recombinant Interferon- α 2c
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human IFN- α subtypes 2a, 2b and 2c
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of human IFN- α subtypes 2a, 2b and 2c in solution e.g. cell culture supernatants using ELISA and for enumeration of IFN- α subtype 2a, 2b and 2c producing cells using ELISpot. MT1 is recommended as coating mAb in combination with detection mAb MT2 (product code 3423-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.