

Monoclonal antibody to Human Interleukin-10

Antibody:	9D7
Product code:	3430-3-1000
Quantity:	1000 µg
Immunogen:	Recombinant human IL-10
Isotype:	Rat IgG1
Specificity:	Native and recombinant human IL-10
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of human IL-10 in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-10 producing cells using ELISpot. 9D7 is recommended as coating mAb in combination with detection mAb 12G8 (product code 3430-6). Immunostaining of human IL-10.
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-04-07



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com