

# Monoclonal antibody to Mouse Interleukin-10

<b>Antibody:</b>	16E3
<b>Product code:</b>	3431-3-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Recombinant mouse IL-10
<b>Isotype:</b>	Rat IgG2b
<b>Specificity:</b>	Native and recombinant mouse IL-10
<b>Concentration:</b>	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For enumeration of IL-10 producing cells using ELISpot. 16E3 is recommended as coating mAb in combination with detection mAb 2A5 (product code 3431-6).
<b>Protocols:</b>	Please contact Mabtech for ELISpot protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on 2020-04-07



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)