

Monoclonal antibody to Mouse Interleukin-10

BIOTINYLATED

Antibody:	16E3
Product code:	3431-4-1000
Quantity:	1000 µg
Immunogen:	Recombinant mouse IL-10
Isotype:	Rat IgG2b
Specificity:	Native and recombinant mouse IL-10
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Applications:	For quantification of IL-10 in solution, e.g. cell culture supernatants, using ELISA. 2A5 (product code 3431-5) is recommended for coating and 16E3 for detection.
Protocols:	Please contact Mabtech for ELISA protocols.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Intertek



Intertek

MABTECH AB
Box 1233
SE-131 28 Nacka Strand
Sweden
Tel +46 8 716 27 00
Fax +46 8 716 27 01
mabtech@mabtech.com

MABTECH Inc.
M.E.B. Suite 220, 3814 West Street
OH 45227 Cincinnati
USA
Tel: +1 513 871 4500
Fax: +1 513 871 7353
mabtech.usa@mabtech.com
Toll Free: +1 866 ELI-SPOT

MABTECH Australia Pty Ltd
Australia
Tel: +61 3 9470 4704
Fax: +61 3 8678 3216
mabtech.au@mabtech.com

MABTECH AB
Büro Deutschland
Germany
Tel: +49 40 4135 7935
Fax: +49 40 4135 7945
mabtech.de@mabtech.com

MABTECH AB
Bureau de liaison
France
Tel: +33 4 92 38 80 70
Fax: +33 4 92 38 80 71
mabtech.fr@mabtech.com