

Monoclonal antibody to Human Interleukin-3

Antibody:	IL3A
Product code:	3435-3-1000
Quantity:	1000 µg
Immunogen:	Recombinant human IL-3
Isotype:	Rat IgG1
Specificity:	Native and recombinant human IL-3
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IL-3 in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-3 producing cells using ELISpot. IL3A is recommended as coating mAb in combination with detection mAb IL3B (product code 3435-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com