

# Monoclonal antibody to Human Interleukin-2

<b>Antibody:</b>	17H12
<b>Product code:</b>	3440-0N-500
<b>Quantity:</b>	500 µg
<b>Immunogen:</b>	Recombinant human IL-2
<b>Isotype:</b>	Rat IgG2a
<b>Specificity:</b>	Native and recombinant human IL-2. Cross-reacts with native IL-2 from rhesus and cynomolgus macaques. Cross-reactivity has been tested in ELISA and/or ELISpot. Please contact Mabtech for information on reactivity with other monkey species.
<b>Concentration:</b>	1 mg/ml
<b>Buffer:</b>	Supplied in sterile filtered (0.2 µm) PBS. Contains no preservatives. The endotoxin content in the product is below 10 EU/mg.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
<b>Applications:</b>	Neutralization of human IL-2 bioactivity. It is recommended to establish the optimal antibody concentration for the assay system used.
<b>References:</b>	Jailwala P, et al. Apoptosis of CD4+ CD25(high) T cells in type 1 diabetes may be partially mediated by IL-2 deprivation. PLoS One. 4: 8, 2009.

# MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Updated on 2020-03-13

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)