

Monoclonal antibody to Mouse Interleukin-2

Antibody:	1A12
Product code:	3441-3N-500
Quantity:	500 µg
Immunogen:	Recombinant mouse IL-2
Isotype:	Rat IgG2a
Specificity:	Native and recombinant mouse IL-2
Concentration:	1 mg/ml
Buffer:	Supplied in sterile filtered (0.2 µm) PBS. Contains no preservatives. The endotoxin content in the product is below 10 EU/mg.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Applications:	Neutralization of mouse IL-2 bioactivity. It is recommended to establish the optimal antibody concentration for the assay system used. The antibody is also suitable for immunoassays.
References:	Jones GW, et al. Naive and activated T cells display differential responsiveness to TL1A that affects Th17 generation, maintenance, and proliferation. <i>FASEB J.</i> 25: 1, 2011.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-03-13



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com