

Monoclonal Antibody to Human Interleukin-12/-23 (p40)

BIOTINYLATED

Antibody:	MT618
Product code:	3450-6-250
Quantity:	250 µg
Immunogen:	Recombinant human IL-23
Isotype:	Mouse IgG1
Specificity:	Native and recombinant p40 subunit of human IL-12 and IL-23. The antibody recognises p40 monomer, p40 homodimer, p70 heterodimer (p40+p35) and IL-23 heterodimer (p40+p19). The antibody cross-reacts with native p40 from rhesus and cynomolgus macaques and recombinant p40 from rhesus macaques. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IL-12/-23 (p40) in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-12/-23 (p40) producing cells using ELISpot. MT618 is recommended as detection mAb in combination with coating mAbs MT86/221 (product code 3450-3).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Updated on, 2020-03-26

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com