

# Monoclonal antibody to Mouse Interleukin-12/-23 (p40)

<b>Antibody:</b>	C17.8
<b>Product code:</b>	3456-3-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Recombinant mouse IL-12
<b>Isotype:</b>	Rat IgG2a
<b>Specificity:</b>	Native and recombinant p40 subunit of mouse IL-12 and IL-23. The antibody recognises p40 monomer, p40 homodimer, p70 heterodimer (p40+p35) and IL-23 (p40+p19).
<b>Concentration:</b>	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of IL-12/-23 (p40) in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-12/-23 (p40) producing cells using ELISpot. C17.8 is recommended as coating mAb in combination with detection mAb IL12-IV (product code 3456-6) for detection of p70.
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on 2020-04-07



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)