

Human Interleukin-23 (p19)
Monoclonal antibody
FOR NEUTRALIZATION

General


Antibody:	MT155
Product code:	3457-5N-500
Quantity:	500 µg
Immunogen:	Recombinant human IL-23
Isotype:	Mouse IgG2b
Specificity:	Native and recombinant human IL-23. The antibody recognises the p19 subunit of the IL-23 heterodimer (p40+p19). The antibody cross-reacts with native p19 from cynomolgus macaques. Cross-reactivity has been tested in ELISA or ELISpot. Please contact Mabtech for information on reactivity with other monkey species.
Concentration:	1 mg/ml
Buffer:	Supplied in sterile filtered (0.2 µm) PBS. Contains no preservatives. The endotoxin content in the product is below 10 EU/mg.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Applications:	Neutralization of human IL-23 bioactivity. It is recommended to establish the optimal antibody concentration for the assay system used. The antibody is also suitable for immunoassays.

Neutralization

Neutralization assay: Neutralization of human IL-23 induced IFN- γ production (see graph). To determine the neutralizing capacity and ND₅₀, mixtures of various concentrations of mAb MT155, human recombinant IL-23 (10 ng/ml) and recombinant human IL-18 (50 ng/ml) were added to human PBMC. Following overnight incubation, cell supernatants were quantified in IFN- γ ELISA.

Neutralization dose 50 (ND₅₀): MT155 neutralized the effect of IL-23 at a ND₅₀ value of 0.2-0.8 μ g/ml when 10 ng/ml recombinant human IL-23 was present.

ND₅₀ is defined as the concentration of antibody that results in 50% inhibition of the biological response under the experimental conditions specified.



IL-23 (10 ng/ml)			+	+	+	+	+	+
IL-18 (50 ng/ml)		+		+	+	+	+	+
MT155 (μ g/ml)					0.1	1	10	100

Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.



The products are for research use only.

MABTECH shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages there from.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com